	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.0
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2022.11.01

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多宁/DuoNing

动物细胞高性能培养基 High-Performance Culture Medium for Animal-Cells

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Maxfeed™ 403A1 and Maxfeed™ 403B1 medium are no animal-derived components and no protein composition, chemical defined high-efficiency feed medium. They are suitable for fed-batch culture of Chinese hamster ovary cells (CHO). These products do not contain growth factors, hypoxanthine, thymine and L- glutamine, and are suitable for the culture of CHO-K1, CHO-S, DG44, CHOZN and other different cell lines using GS system and dhfr screening system.

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Maxfeed™ 403A1 should be used together with Maxfeed™ 403 B1 and added simultaneously.

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When using, it is necessary to combine the characteristics of the cell line, and use a reasonable feeding strategy can significantly improve the cell growth and increase the production of antibody expression.

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① Maxfeed™ 403A1: Does not contain glucose, the preparation process can be adjusted according to the actual process, recommended to add 75g / L of glucose. The pH of liquid Maxfeed™ 403A1 after preparation was neutral.

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.0
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2022.11.01

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Maxfeed™ 403B1 Medium: does not contain glucose, the pH of liquid Maxfeed™ 403B1 medium is alkaline.

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② For the first time, it is suggested to start feeding in the middle and late period of cell exponential growth, usually on the third day when the cell density reaches $2.0\sim 3.0\times 10^6$ cells/mL. The recommended feeding scheme is as follows, and can be adjusted according to the cell lines used:

补料时间 Feeding time	Maxfeed™ 403A1+ Maxfeed™ 403B1
Day 3	3%+0.3%
Day 5	3%+0.3%
Day 7	5%+0.5%
Day 9	5%+0.5%
Day 11	3%+0.3%

注意事项 matters need attention:

- 培养 GS 筛选系统工程细胞株时, 根据需要添加次黄嘌呤和胸腺嘧啶。
- When cultivating the engineering cell of GS screening system, hypoxanthine and thymine were added as needed.
- 培养 dhfr 筛选系统工程细胞株时, 根据需要添加 3~8mM L-谷氨酰胺。
- When cultivating dhfr screening system engineering cell lines, 3~8mM L- glutamine was added as needed.

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Suitable for Maxfeed™ 403A1 produced by Wuxi Duoning pin mill process (taking 1L as an example)

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1.Measure 80% of the final volume ultrapure water (20~30 °C);

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3.Add 96.18g of Maxfeed™ 403A1 powder, stir for 30 min, and dissolve completely. At this time, the pH is about 6.0;

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4. Add 75g/L glucose and continue stirring for 10min;

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Remarks: Maxfeed™ 403A1 powder will be turbid when it is initially dissolved, and the solution will gradually clarify with the increase of pH.

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6. Constant volume, continue to stir for 5 ~ 10 min;

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7. Sterilizing and filtering with a 0.22μm filter.

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Suitable for Maxfeed™ 403B1 medium produced by Wuxi Duoning pin mill process (taking 1L as an example)

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1. Measure 70% of the final volume ultrapure water (20~30 °C);

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3. Add 82.7g of Maxfeed™ 403B1 medium, stir for 30min, and dissolve completely;

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Shanghai production base, Maxfeed™ 403A1, powder packaging: store in dark at 2~8 °C, and the validity period is 24 months.

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Wuxi production base, Maxfeed™ 403A1, powder packaging: store in dark at 2~8 °C, and the retest date is 24 months.

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Shanghai production base, Maxfeed™ 403B1 medium, powder packaging: store in dark at 2~8 °C, and the validity period is 24 months.

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Wuxi production base, Maxfeed™ 403B1 medium, powder packaging: store in dark at 2~8 °C, and the retest date is 24 months.

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地址 Address: 上海市松江区民强路 1525 弄 30 幢 Building 30, Lane 1525, Minqiang Road, Songjiang District, Shanghai

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