	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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多宁/DuoNing

动物细胞高性能培养基 High-Performance Culture Medium for Animal-Cells

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DN feed 3 & DN feed B3 are chemical defined (CD) highly efficient feed medium with no animal derived components, no protein, which are suitable for use in fed-batch cultures using Chinese hamster ovary cells (CHO) in the development and production of therapeutic protein products. This product does not contain growth factors, hypoxanthine, thymidine or L-glutamine. It is suitable for the culture of different cell lines such as CHO-K1, CHO-S, DG44, CHOZN using GS system and DHFR screening system.

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DN feed 3 should be used together with DN feed B3 and added at the same time.

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DN feed 3 & DN feed 4 are highly efficient supplementation formula designed to improve the growth and production performance of CHO cells. They are rich in amino acids and vitamins to meet the nutritional requirements of high-density cell culture and high product expression. Combined with the cell line characteristics, a reasonable supplementation strategy can significantly improve the cell growth and antibody expression performance.

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The first feeding is recommended in the middle to late stage of cell index growth, usually on day 3 when the cell density reaches $2.0-3.0 \times 10^6$ cells/mL or more, and the recommended feeding scheme is as follows, and the optimized ratio can be customized according to the cell line used:

补料时间 Feeding time	DN feed 3+DN feed B3
Day 3	3%+0.3%
Day 5	3%+0.3%
Day 7	3%+0.3%
Day 9	3%+0.3%
Day 11	3%+0.3%

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DN feed 3 Liquid packaging: 2°C to 8°C, protect from light; validity period: 12 months.

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DN feed B3 Liquid packaging: 2°C to 8°C, protect from light; validity period: 12 months.

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