	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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多宁/DuoNing

动物细胞高性能培养基 High-Performance Culture Medium for Animal-Cells

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Transpro feed 1 medium is a universal transient transfection feed medium that can be used in conjunction with Transpro CD 01 medium. This product can be used for transient transfection culture of both HEK293 and CHO cells. Transpro feed 1 is suitable for transient transfection and expression culture of antibodies, recombinant proteins, and viruses during the research and development process using HEK293 series cells such as HEK 293, Expi293F, 293F, 293E, and CHO series cells such as ExpiCHOS and CHOS. This product is a fully chemically defined medium, without animal derived ingredients, protein derived ingredients, animal or plant derived protein hydrolysates, or growth factors.

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Transpro feed 1 contains a certain concentration of transient enhancer, which is suitable for single feeding. The recommended feeding volume is 3-10% of the initial cultivation volume (recommended 5%). This product is not recommended to be supplied for many times, and it is not necessary to add transient transfection enhancer.

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Cell transient transfection (using HEK293 transient protein as an example)

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The day before transfection need to seed cells at 2.0×10^6 cells/mL, the cell density can reach 4.0×10^6 cells/mL on the second day;

	记录编号 FN: SMP-RD-A0002-R06	版本号 Ver: 2.1
	归档 Filing: filed immediately, permanently kept	生效日期 Ed: 2023.01.12

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After cell counting on the second day of culture, the cell viability was more than 95%, and the living cell density was $\geq 4.0 \times 10^6$ cells / mL, can be used directly; If the cell density is lower than 4.0×10^6 cells / mL, the cells can be collected by centrifugation (800 rpm, 5 min), and the cells can be resuspended at 4.0×10^6 cells/mL density in Transpro CD 01;

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The mixture of DNA and PEI was prepared according to the optimized transient process;

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Add the mixed solution to the medium for culture;

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After 18 hours of culture, it is recommended to supplement the supplemented medium Transpro feed 1 (the concentration is recommended to be 3-5% of the initial culture volume), or the combined supplemented medium DN feed B2 (the concentration is recommended to be 0.3-0.5% of the initial culture volume), which can further improve the density of living cells and protein production.

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Detecting residual sugar in the process of cell culture medium, and controlling the concentration of residual sugar at 3-6g/L;

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Culture until 7 days, or the viability is less than 60%, and end the culture.

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Suitable for powder packaging (taking 1L as an example).

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Slowly add 200.879 g of Transpro feed 1 medium powder, stir for 30 min and dissolve completely;

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Constant volume, stirring for 5 ~ 10 min, and detecting osmotic pressure (osmotic pressure range is 2000~2400 mOsm/kg);

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Sterilization and filtration by 0.22 μ m filter.

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Wuxi production base, Transpro feed 1 medium powder packaging: 2 $^{\circ}$ C to 8 $^{\circ}$ C, protect from light; retest date: 24 months.

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